

#9463 Gel Beads Hot/Cold Pack



Regular Pricing

| | | | | | |
|--------|--------|--------|--------|--------|----|
| 100 | 250 | 500 | 1000 | 2500 | |
| \$4.12 | \$3.58 | \$3.11 | \$2.70 | \$2.35 | 5C |

Highlights

- Therapeutic Gel Pack Applies Heat Or Cold To Sore Muscles
- Microwave And Freezer Safe
- Reusable And Non-Toxic
- Instructions Printed On Reverse Side

Packaging Details

| Quantity per Box | Box Weight | Box Length | Box Width | Box Height |
|------------------|------------|------------|-----------|------------|
| 50 | 23 lbs. | 13" | 13" | 11" |

Dimensions are in inches. Weight is in pounds.

Description

- COLORS AVAILABLE: Blue, Clear, Green, Purple, Red, Dark Blue, Pink or White.
- IMPRINT COLORS: Standard Pad-Print Colors
- APPROXIMATE SIZE: 7 1/4" W x 4 1/4" H
- IMPRINT AREA AND METHOD: 3" W x 2" H
- SET UP CHARGE: \$40.00(G), \$25.00(G) on re-orders.
- MULTI-COLOR IMPRINT: Not Available
- PACKAGING: Poly Bag
- PLEASE NOTE: After freezing, stays cold for up to 3 hours. If dropped, the item will not shatter or break. After microwaving 30 seconds, the item will stay warm for 90 minutes. Condensation may potentially cause bubbling until cool. CAUTION: Product will be HOT after microwaving.